



Efficient Wound Healing Without Disturbance!

SERIDERM[®] BORDER – For Smoother and Swifter Healing



Polyurethane Layer

- Protective barrier against bacteria and contaminants
- Prevents wound infections
- Maintains a moist environment around the wound for optimal healing
- Breathable and waterproof

Silicone Adhesive

- Adheres to the skin around the wound
- Forms a secure seal without disturbing the healing process
- Painless application and removal
- Reduces periwound complications

Cotton Gauze with 100% Silk Protein in Wound Contact Layer

- Soft, absorbent layer
- Wicks away excess exudate for optimal moisture management
- Improvement in inflammatory markers
- Supports tissue regeneration for expedited healing

Indications

SERIDERM[®]BORDER is the ideal wound dressing for a wide range of chronic and acute wounds that are infected or high-risk wounds, including:


- Diabetic ulcers
- Surgical wounds
- Pressure ulcers (bed sores)
- Wound preparation for grafting
- Donor Sites

Also Available As:



SERIDERM[®] Ag BORDER:
Superior Healing of SERIDERM[®]BORDER
Meets the Antimicrobial Action of
Silver Ions!

**For Infected or High-Risk
Wounds**



ORDER NOW



+91-8468847349



info@serigenmed.com

SERIDERM [®] BORDER	Size (in cm)	Item Code
Square	7.5x7.5	SDB-7575
Square	10x10	SDB-1010
Square	15x15	SDB-1515
Square	20x20	SDB-2020
Square	23x23	SDB-2323

SERIDERM [®] Ag BORDER	Size (in cm)	Item Code
Square	7.5x7.5	SAB-7575
Square	10x10	SAB-1010
Square	15x15	SAB-1515
Square	20x20	SAB-2020
Square	23x23	SAB-2323



SERIGEN MEDIPRODUCTS PVT. LTD.

39, Electronic Co Operative Estate Parvati Paytha,

Pune, Maharashtra 411009, India

www.serigenmed.com

Version - 05/MAR/2025