





The Power of Two to Combat Wound Infections!

Seriderm[®]Ag – A Complete Wound Healing Solution For Infected or High-Risk Wounds

1. Healing power of Silk Protein



2. Antimicrobial action of Silver Ions



 **SERIGEN**




100% silk protein in wound contact layer

Sustained release of silver ions for up to 7 days



Seriderm[®] + Silver = Seriderm[®] Ag
A complete wound-healing solution!




Minimize Wound Infections with Silver!


3-Way Antimicrobial Action of Silver Against Infections



Disturbing the permeability of the bacterial cell membrane



Interaction with biomolecules inside the microbial cells



Interruption of the respiratory chain



Seriderm[®] Ag is a non-adherent, absorbent wound dressing designed using **98% silk protein and silver ions**.

Seriderm[®] Ag offers **accelerated wound healing and combats infections** in a wide range of infected or high-risk wounds, including diabetic foot ulcers, bed sores, and second-degree burns.

✓ Faster Healing Powered By Silk

The novel silk-based composition of Seriderm[®] Ag modulates cellular pathways resulting in better vasculature and granulation tissue formation for expedited wound healing and reduction of inflammatory markers in wounds.

✓ Moist Wound Healing Environment

Seriderm[®] Ag can absorb exudate >900% of its weight and maintains optimal moisture at the wound environment for accelerated healing, with a moisture vapour transmission rate >1000 g/m²/24h.

✓ Sustained Silver Release

Seriderm[®] Ag is designed for the sustained release of silver ions over 7 days, extending antimicrobial protection and reducing the need for frequent dressing changes.

✓ Broad Spectrum Antimicrobial Protection

Seriderm[®] Ag exhibits broad-spectrum antimicrobial activity against a wide range of wound-related microorganisms, including multi-drug resistant organisms.



Indications

Seriderm[®] Ag is the ideal wound dressing for a wide range of chronic and acute wounds that are infected or high-risk wounds, including:

- Diabetic ulcers
- Surgical wounds
- Pressure ulcers (bed sores)
- Venous ulcers
- Wound preparation for grafting
- Donor Sites
- Second-degree burns

Manufactured in an ISO
13485-certified facility



ORDER NOW

+91-8468847349

info@serigenmed.com

Type	Size (in cm)	Item Code
Square	5x5	SD-AG-0505
Square	10x10	SD-AG-1010
Square	20x20	SD-AG-2020
Square	23x23	SD-AG-2323
Square	30x30	SD-AG-3030



SERIGEN MEDIPRODUCTS PVT. LTD.

39, Electronic Co Operative Estate Parvati Paytha,
Pune, Maharashtra 411009, India

www.serigenmed.com

ver 2.0, July 2026